



SPECIMEN COLLECTION CHART

Tube Type	Specimen	Tests
<p>Red</p>  <p>No Additive</p>	Serum	<p><i>Histoplasma</i> Ag, EIA <i>Blastomyces</i> Ag, EIA <i>Coccidioides</i> Ag, EIA <i>Histoplasma</i> Canine IgG Ab, EIA <i>Histoplasma</i> Feline IgG Ab, EIA <i>Coccidioides</i> Canine IgG Ab, EIA <i>Blastomyces</i> Canine IgG Ab, EIA</p> <p><i>Aspergillus</i> Galactomannan <i>Aspergillus</i> Ab, ID (1→3)-β-D-Glucan <i>Cryptococcal</i> Ag <i>Coccidioides</i> Ab, ID Itraconazole</p>
<p>Gold****</p>  <p>Clot Activator with Gel</p>	Serum	<p><i>Histoplasma</i> Ag, EIA <i>Blastomyces</i> Ag, EIA <i>Coccidioides</i> Ag, EIA <i>Histoplasma</i> Canine IgG Ab, EIA <i>Histoplasma</i> Feline IgG Ab, EIA <i>Coccidioides</i> Canine IgG Ab, EIA <i>Blastomyces</i> Canine IgG Ab, EIA</p> <p><i>Aspergillus</i> Galactomannan <i>Aspergillus</i> Ab, ID (1→3)-β-D-Glucan <i>Cryptococcal</i> Ag <i>Coccidioides</i> Ab, ID Itraconazole</p>
<p>Tiger Top****</p>  <p>Clot Activator with Gel</p>	Serum	<p><i>Histoplasma</i> Ag, EIA <i>Blastomyces</i> Ag, EIA <i>Coccidioides</i> Ag, EIA <i>Histoplasma</i> Canine IgG Ab, EIA <i>Histoplasma</i> Feline IgG Ab, EIA <i>Coccidioides</i> Canine IgG Ab, EIA <i>Blastomyces</i> Canine IgG Ab, EIA</p> <p><i>Aspergillus</i> Galactomannan <i>Aspergillus</i> Ab, ID (1→3)-β-D-Glucan <i>Cryptococcal</i> Ag <i>Coccidioides</i> Ab, ID Itraconazole</p>
<p>Purple</p>  <p>K2 EDTA</p>	Plasma	<p><i>Histoplasma</i> Ag, EIA <i>Blastomyces</i> Ag, EIA <i>Coccidioides</i> Ag, EIA</p> <p><i>Aspergillus</i> Galactomannan</p>

* Label specimens with animal's name and owner/facility name.

** Ship to arrive Monday – Saturday.

*** *Histoplasma*, *Blastomyces* and *Coccidioides* Antigen specimens may be shipped Room Temperature. Specimens for all other testing should be refrigerated or frozen and shipped on cold packs.

**** Specimens can be sent in a serum separator tube or transferred to a plastic screw top tube or red top tube for transit.

HEADQUARTERS

4705 Decatur Blvd. | Indianapolis, Indiana 46241, USA

888-841-8387



SPECIMEN COLLECTION CHART

Tube Type	Specimen	Tests
<p>Dark Green</p>  <p>Sodium Heparin</p>	Plasma	<p><i>Histoplasma</i> Ag, EIA <i>Blastomyces</i> Ag, EIA <i>Coccidioides</i> Ag, EIA</p> <p><i>Aspergillus</i> Galactomannan</p>
<p>Light Green</p>  <p>Lithium Heparin</p>	Plasma	<p><i>Histoplasma</i> Ag, EIA <i>Blastomyces</i> Ag, EIA <i>Coccidioides</i> Ag, EIA</p> <p><i>Aspergillus</i> Galactomannan</p>
<p>Light Blue</p>  <p>Sodium Citrate</p>	Plasma	<p><i>Histoplasma</i> Ag, EIA <i>Blastomyces</i> Ag, EIA <i>Coccidioides</i> Ag, EIA</p> <p><i>Aspergillus</i> Galactomannan</p>
<p>Urine tube/cup</p> 	Urine	<p><i>Histoplasma</i> Ag, EIA <i>Blastomyces</i> Ag, EIA <i>Coccidioides</i> Ag, EIA</p> <p><i>Aspergillus</i> Galactomannan</p>
<p>Dark Green</p> 	CSF, BAL	<p><i>Histoplasma</i> Ag, EIA <i>Blastomyces</i> Ag, EIA <i>Coccidioides</i> Ag, EIA <i>Histoplasma</i> Canine IgG Ab, EIA <i>Histoplasma</i> Feline IgG Ab, EIA</p> <p><i>Coccidioides</i> Canine IgG Ab, EIA <i>Blastomyces</i> Canine IgG Ab, EIA <i>Aspergillus</i> Galactomannan (1→3)-β-D-Glucan Cryptococcal Ag</p>

HEADQUARTERS

4705 Decatur Blvd. | Indianapolis, Indiana 46241, USA

888-841-8387